

**Doctor recommended** Laser Hair Growth Treatment for Safe and Effective Home Use

## What is Theradome®?

- Wearable laser phototherapy (LPT) hair growth treatment
- Reduces hair loss
- Thickens existing hair follicles by reversing miniaturization
- Grows new hair follicles
- The only device with the most advanced VL680 laser diodes designed to last 18,000 treatments

## Who is it for?




- Those affected by androgenetic alopecia who want their hair back and to stop hair loss
- People that want to treat hair loss in the comfort of their own home with zero side effects.
- FDA-cleared for the safe and effective treatment of female and male pattern baldness

## How does it work?

- Only Theradome's clinical strength patented VL680 lasers, provide the optimal dose in minimal time
- Comfortable helmet design that covers more of the scalp than *any other device*
- Targets the bulge and the bulb at the base of hair follicles for maximum mitochondrial stimulation
  - Bulge maintains hair growth cycle keeping it in the active, healthy anagen phase
  - Bulb produces Keratin for the generation and speed of hair growth


## Why Theradome®?

- 98% of patients saw improved hair growth and reduced hair loss using Theradome
- Theradome® is more powerful because of the type of lasers used that are made in the USA
- Theradome® lasers reaching hair follicle base unlike EEL and LED that are sourced in Asia
- Use for 20 minutes, 2 x per week
- Heat management system. Uninterrupted laser dosage delivered without overheating
- Multi-award winning technology



# Theradome® LH80 PRO hair restoration device

Technologically Advanced Laser Hair Growth Treatment for safe and effective home use. Clinical trials have shown that hair loss can be reversed with the Theradome® LH80 PRO Laser Helmet. Result can be seen in as soon as 4 weeks although this may vary from person to person.



See results with the Theradome® LH80 PRO in as soon as 4 weeks!

- ✓ Flexible + convenient
- ✓ Simple to use
- ✓ Hands-free
- ✓ Clinically tested
- ✓ At home use
- ✓ Light weight

theradome